

Guidance Sheet 12:

Lung Function Sub Study



For Sites Taking Part in the Lung Function Sub Studies

Introduction

The Tremoflo N-100 Airwave Oscillometry System is a portable medical device intended to monitor lung function in newborns and infants. The device measures airway resistance, reactance and other lung function parameters.

Although the device is portable, please only use it in the hospital and is not intended for home use. The device can be used for research purposes only and is not to be used for diagnostics.

Equipment

You will be provided with the following equipment:

- Tremoflo N-100 Airwave Oscillometry System
- Laptop with the Tremoflo software installed
- Masks (3 different sizes) (single use)
- Filters (single use)

Set Up

The laptop will arrive with the necessary software installed and activated to use the Tremflo. The following equipment will be provided:

Guidance Sheet 12:

Lung Function Sub Study



In addition, a power supply and an ethernet cable will be supplied. A calibration test load and test load adapter will also be supplied.

The Tremflo will need calibrating, to calibrate the machine please contact the BALLOON Trial Manager (kotechasj@cardiff.ac.uk) to set up a Team's call at a convenient time with Thorasys, the Tremflo manufacturer. This call needs to take place in the afternoon (as they are based in Canada) and please allow 2 hours for the calibration, although it may not take that long.


Taking a Measurement

To power up the Tremoflo system, briefly press the power button on the Tremoflo cradle. The power button will light up, briefly alternating between green, yellow and red before settling to a steady yellow colour.

The Tremoflo has two power indicator lights. On the cradle, a power and status indicator light around the power button can light up in green, yellow or red, each colour has a significance as shown in the below table. The second power indicator, located on the top of the handheld unit, lights up in blue when power is supplied to the handheld unit and does not change colour with the system status.

Guidance Sheet 12:

Lung Function Sub Study




Colour	Status	Details
Unlit/Gray	Off	The tremoflo is powered down.
Flashing	Starting	The tremoflo firmware is booting.
Green	Ready	The tremoflo is powered up, communicating with the software, fully initialized and ready for use.
Yellow	Standby	The tremoflo is powered up and either waiting to connect to the software or being initialized.
Red	Error	The tremoflo is powered up, but a critical error has occurred and the system is not ready for use. Refer to Appendix A for possible causes and solutions.

With the computer running start the Tremoflo software by either clicking the Tremoflo icon in the task bar or by double clicking the Tremoflo icon of the desktop.

Enter your username and password and then click login.

Guidance Sheet 12:


Lung Function Sub Study



Next the software shows a Welcome dialog.

New Patient: Click to add and test a new patient.

Select patient: To select an existing patient from the database to retest or review collected data. Please see note below *



To add a new patient select **New Patient** from the **Welcome** dialog or from the Tremoflo toolbar on the **Patient** tab to display the **New Patient** dialog. Enter the patient information.

Guidance Sheet 12:

Lung Function Sub Study



Field	Information Entered	Format
First Name	Patient's first name	Alphanumeric Text
Last Name	Patient's last name	Alphanumeric Text
Subject Code	Custom Code	Alphanumeric Text
Sex	Male/Female	From drop-down list
Date of Birth	Date Of Birth	As per regional settings
Weight	Body Weight	See Table 4.3
Height	Body Height	See Table 4.3
Ethnicity	Patient's Ethnicity	From drop-down list

Please use centimetres for length/height and kilograms for weight.

If you wish to choose a patient who is already in the database please choose select **Patient**. *Please enter a participant as new patient when you see them at baseline, and another new patient at 6 months corrected age and another new patient at 12 months corrected age. When you see the participant at baseline please use the PID with a 0 at the end . When you see the participant at 6 months of age please use the PID with a 6 at the end. When you see the participant at 12 months of age please use the PID with a 12 at the end.

Guidance Sheet 12:

Lung Function Sub Study



tf Patient List

All Patients

	Jane Doe Birthday: 19/10/1988 Age: 27.0 yrs Last Tested: Never	Sex: Female Weight: 58 kg Height: 155 cm
	John Doe Birthday: 19/10/1977 Age: 38.0 yrs Last Tested: 19/10/2015 12:36:25 PM	Sex: Male Weight: 75 kg Height: 165 cm

Test Patient Review Data Cancel

Guidance Sheet 12:

Lung Function Sub Study



tremoFlo 1 | Patient 2 | Test 3 | Results SUPPORT MODE K, D

New Patient Select Patient Edit Patient New Test View Charts View Table Patient Details

Date of Birth: Oct. 14, 1991
Age: 31.1 yrs
Weight: 64 kg
Height: 165 cm
Last Tested: Oct. 14, 2022

14/10/22 Time

14/10/22 Time

14/10/22 Time

14/10/22 Time

D K

Identification

TFID 194

Added 10/14/2022 8:20:16 PM

Subject Code

Family Name K

Given Name D

Patient Details

Sex Female

Age 31.1 yrs

Date of Birth 10/14/1991

Weight 64 kg

Height 165 cm

Ethnicity Unspecified

Smoker


Pack Years 0

OK Cancel

When choosing **“Test a patient”** the test type dialog will appear with a choice of two test types:

Guidance Sheet 12:

Lung Function Sub Study



Please select **Standard test**.

Once you select **Standard test** type the **Test properties** dialog is displayed. Chose the **Reference data set** (You will choose N-100 template) and fill in the other fields as required.

Guidance Sheet 12:

Lung Function Sub Study



New Test

Title: Test 1

Template: Airwave Oscillometry

Contains waveforms:
Airwave Oscillometry 5-37 Hz
SoftPulse Oscillometry 5-25 Hz

Reference Dataset: Oostveen E. et al. ERJ Express 2013

Consult User Manual for further information on the use of reference data.

Operator: admin

Attending Physician:

Prior Medication:

Operator Comments:

Status: Open

Skip this dialog

OK Cancel


Closing the **Test properties** dialog completes the start-up sequence of the Tremoflo software. At this point, the device will connect to the software and be initialised. Once the system initialisation is complete the screen will display the main window of the Tremoflo software in the **Test Tab**. At this point, you are ready to take your first

Guidance Sheet 12:

Lung Function Sub Study



measurement. The machine may ask you to calibrate it, please use the 90 load to do this and attach it to the end of the tubing using the load adapter.



Please now prepare the baby for testing. The baby should be asleep and if possible swaddled to prevent movement.

1. Make sure the Tremoflo is showing the test screen.
2. Use gloves when handling the Tremoflo
3. Remove a new and packaged filter and mask from its wrapping and insert into the machine as show in the diagram below.
4. Position the baby as shown below

Guidance Sheet 12:

Lung Function Sub Study



(link to video on website)

The video at <https://www.youtube.com/watch?v=umMhsT-6ke4> shows the machine in action in a child.

Once the baby is breathing quietly, please complete at least **three** satisfactory measurements of 25-30 seconds each by pushing the button on the handheld unit to start the measurement. The device will vibrate and you will hear a quiet noise. The Tremoflo automatically ceases measuring when the measurement is complete. Measurement data are automatically stored without the need for user intervention.

When the test is complete, the test should be marked as closed in the Tremoflo software.

Please dispose of the mask and filter following your Trust/Health Board guidance. It is important to clean the exterior of the machine using a mild detergent of a neutral pH with warm water and a soft cloth. Take care not to press against or otherwise damage the mesh screen inside.

Downloading the data

Please download the data after each patient by following these steps and transfer securely to the BALLOON team.

Guidance Sheet 12:


Lung Function Sub Study



Go to the Tremoflo menu in the software (top left corner).

Select the “**Export Data**” option.

On the **Data Export Wizard**, select the needed export scenarios and the export format.



Click on the “**Next**” button.

Select the patients containing the data that need to be exported. Click “**Next**”.

Guidance Sheet 12:

Lung Function Sub Study



Data Export

Step 1: Scenario Selection
Step 2: Patient Selection
Step 3: Patient Test Selection
Step 4: Export Timeframe
Step 5: Export Destination
Step 6: Export Process

Step 2: Patient Selection

Identification	Sex	Birthdate	Weight	Height
<input checked="" type="checkbox"/> 123 - S D	Male	10/14/1996	79 kg	172 cm
<input type="checkbox"/> 123 - P J	Female	9/5/1995	65 kg	180 cm
<input type="checkbox"/> D K	Female	10/14/1991	44 kg	165 cm
<input type="checkbox"/> F M	Female	9/21/2016	30 kg	140 cm
<input type="checkbox"/> Flou M	Female	11/18/2002	59 kg	163 cm
<input type="checkbox"/> François Morin	Male	12/21/1976	77.112 kg	176 cm
<input type="checkbox"/> T P	Male	11/7/2000	75 kg	180 cm
<input checked="" type="checkbox"/> Test Patient	Male	9/12/2003	78 kg	180 cm
<input type="checkbox"/> N T	Male	9/22/1995	85 kg	176 cm
<input type="checkbox"/> Keith Westcott	Male	3/17/1970	84.37 kg	177.8 cm

All
None
Current

< Back Next > Cancel

For each Patient, select the tests containing the data that need to be exported. Click “Next”.

Guidance Sheet 12:

Lung Function Sub Study



Data Export

Step 3: Patient Test Selection

Step 1: Scenario Selection
Step 2: Patient Selection
Step 3: Patient Test Selection
Step 4: Export Timeframe
Step 5: Export Destination
Step 6: Export Process

ID	Title	Timestamp	Measurement	
Test Patient, 19.2 yrs, Male, 78 kg, 180 cm, Unspecified				
<input checked="" type="checkbox"/>	793	Bronchodilator Test: Pre BD	9/12/2022 10:05:22 AM	3
<input checked="" type="checkbox"/>	794	Bronchodilator Test: Post BD	9/12/2022 10:16:04 AM	3
<input checked="" type="checkbox"/>	799	Standard Test	9/19/2022 3:59:22 PM	3
<input type="checkbox"/>	800	Standard Test	9/19/2022 4:05:27 PM	3
<input type="checkbox"/>	801	Standard Test	9/19/2022 4:09:14 PM	1
<input type="checkbox"/>	802	Bronchodilator Test: Pre BD	9/21/2022 9:52:05 AM	3
<input type="checkbox"/>	803	Bronchodilator Test: Post BD	9/21/2022 10:00:52 AM	0
<input type="checkbox"/>	804	Standard Test	9/21/2022 2:42:40 PM	1
<input type="checkbox"/>	805	Standard Test	9/21/2022 2:51:48 PM	0
123 - S D, 26.1 yrs, Male, 79 kg, 172 cm, Unspecified				
<input checked="" type="checkbox"/>	817	Standard Test	10/14/2022 8:12:28 PM	3

All
None
Current

< Back Next > Cancel

Select the timeframe. Click **“Next”**.

Guidance Sheet 12:

Lung Function Sub Study



The screenshot shows a software window titled "Data Export" with a sidebar on the left containing a list of steps: Step 1: Scenario Selection, Step 2: Patient Selection, Step 3: Patient Test Selection, Step 4: Export Timeframe (highlighted with a blue square), Step 5: Export Destination, and Step 6: Export Process. The main area is titled "Step 4: Export Timeframe" and contains two radio button options: "Entire Patient timeframe" (selected) and "Between specific events". Below these are two dropdown menus for "First:" and "Last:". The "First:" dropdown shows "9/12/2022 10:10:10 AM Measurement: AOS 5-37" and the "Last:" dropdown shows "10/14/2022 8:15:10 PM Measurement: AOS 5-37". At the bottom right, there are three buttons: "< Back", "Next >" (highlighted with a blue border), and "Cancel".

Choose the default exportation path or switch it with your own one. Click **“Export”**.

Guidance Sheet 12:

Lung Function Sub Study



Data Export

Step 1: Scenario Selection
Step 2: Patient Selection
Step 3: Patient Test Selection
Step 4: Export Timeframe
Step 5: Export Destination
Step 6: Export Process

Step 5: Export Destination

Export destination folder
C:\Users\Public\Documents\THORASYS\tremoFlo 1.0\ Browse

If File Exists... Prompt

< Back Export Cancel

If the participant is for training only, please also send the following demographics.

- BALLOON ID
- Birthweight
- Current Length
- Current Weight
- Gestation at birth
- CGA
- Date of measurement
- DOB Gender
- On supplementary oxygen

Guidance Sheet 12:

Lung Function Sub Study




To transfer securely

Save the document as a PDF

Visit <https://fastfile.cardiff.ac.uk/>

Fig 1.



Guidance Sheet 12:

Lung Function Sub Study



Fig 2.

Fast File

Home Login

Information about the Sender

If you have been given a "Request Code" then just enter it here and click the button at the bottom of this form.

Request Code:

If you do not have a "Request Code" then please complete the rest of this form:

Your name: (required)

Your organisation: (required)

Your email address: (required)

I now need to send you a confirmation email.
When you get it in a minute or two, click on the link in it.

Please complete these fields

Press 'send' for confirmation.

Fig3.

Fast File

Home Login

Now wait for the email message from the FastFile to arrive and click on the link in it.

You may close this window.

You will be directed to the [main menu](#) in a moment.

Please check your mailbox and follow the link

FastFile User's Guide V1.0 2010-05-26 (Page 1)

Guidance Sheet 12:

Lung Function Sub Study



Fig 4.

Please click on this symbol

Enter BALLOON

Enter BALLOON@cardiff.ac.uk

Guidance Sheet 12:

Lung Function Sub Study



Choose the File(s) you would like to upload

File 1: OASIS Trial - PID 3009.docx
Description: (1.91MB) ✖

File 2:
 No file chosen
Description: ✖

File 3:
 No file chosen
Description: ✖

1.91MB /

15360MB

Choose the file to be upload

Click here.

Troubleshooting

If you have any problems using the Tremoflo system please contact the BALLOON Trial team (BALLOON@cardiff.ac.uk). We will be able to put you in touch with Thorasys' technical support if we are unable to help you.